CA-MI

Kiwi+

Dispositivo per Aerosolterapia

FULL OPTIONAL

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KIWI Plus is an aerosoltherapy device operating with new generation MESH technology. The principle of nebulisation on which it is based, combines, in a single device, the best features of compressor technology (rapid and effective nebulisation of any type of drug) and ultrasonic technology (absolute silence and dual power supply, mains or battery), making it ideal for any type of use on adults and children, both in home-care and in emergency, while travelling or outdoors. All this in addition to construction/operational features developed in the design phase and unique to this type of equipment. The drug containment and nebulisation chamber is made of PC (polycarbonate), a material that ensures high resistance. Over 200 minutes of autonomy with internal rechargeable Lithium-Iones battery and the special design of the nebulisation chamber ensures maximum therapeutic efficacy, reducing the amount of non-spray composition drug to less than 0.1 ml (the lowest declared on the market). The mesh technology in this product has been optimized to the higher level, making possible to achieve particularly high and effective nebulization parameters with **93% breathable fraction and an MMAD of 1.51 µm**

CAUTION

CONSULT THE USER MANUAL BEFORE USING THE DEVICE

ALWAYS FOLLOW YOUR DOCTOR'S INSTRUCTIONS WHEN ADMINISTERING THE DRUG

NEVER DISASSEMBLE THE DEVICE. FOR ANY INTERVENTION CONTACT THE DISTRIBUTOR AND/OR CA-MI TECHNICAL ASSISTANCE

REMOVE THE DEVICE FROM THE CARRYING CASE BEFORE EACH USE

GENERAL WARNINGS

1. When opening the package, check the integrity of the device and the power supply. Pay particular attention to any damage to the plastic parts, which may expose live inner parts, and to breaks and/or peeling of the power supply and relative power cable. **In such cases, do not connect the plug to the power socket. Always carry out these checks before each use.**
2. Before connecting the device, always check that the electrical data indicated on the power supply data label and the type of plug used, match those of the mains electric power supply to which it is to be connected.
3. Follow the safety regulations for electrical equipment and in particular:
 - Exclusively use original Accessories and Spare parts supplied by Manufacturer. Only use original power supplies of the type indicated in this manual.
 - **DO NOT THROW DEVICE IN THE FIRE OR PLACE NEAR FLAMES OR HEAT SOURCES. THE DEVICE MAY EXPLODE.**
 - Never immerse the device and power supply in water;
 - Put the device on a flat and stable surface in order to avoid accidental falls or overheating the device or power supply;
 - Do not use the device in environments where there are anaesthetic mixtures inflammable with air, oxygen or nitrous oxide;
 - Do not use the device in anaesthesia and pulmonary ventilation systems
 - Do not touch the device with wet hands and always prevent the device from coming into contact with liquids;
 - The use of this device by children and/or physically or mentally impaired people always requires the careful supervision of an adult with full mental faculties;
 - Do not leave the device on and/or connected to the power socket when it is not used;

- Do not pull the power supply cable to unplug it, but grasp it with your fingers to extract it from the mains power socket;
 - Store and use the device and all its parts in places protected against the weather and far from any sources of heat. After each use, it is recommended to store the device in its carrying case away from dust and sunlight.
 - In general, it is inadvisable to use single or multiple socket adapters and/or extensions. Should their use be necessary, you must use ones that are in compliance with safety regulations, however, taking care not to exceed the maximum power supply tolerated, which is indicated on the adapters and extensions.
4. For repairs, contact only the Manufacturer technical service department or technical service centre authorised by it and request the use of original spare parts. Failure to comply with the above can compromise the safety of the device
 5. **This device must be used exclusively for its intended use and as described in this manual. Therefore, it must be used as an aerosoltherapy system.** Any use differing from the intended use is considered to be improper and therefore dangerous and neither the Manufacturer nor the Distributor can be held responsible for damages caused by improper, incorrect and/or unreasonable use or by the connection to electrical systems not compliant with current safety standards.
 6. The medical device requires special precautions regarding electromagnetic compatibility and must be installed and used in accordance with the information provided with the accompanying documents. The **KIWI Plus** device must be installed and used away from mobile and portable RF communication devices (mobile phones, transceivers, etc.) that may interfere with it.
 7. Children and dependent people must always use the medical device under the strict supervision of an adult in full possession of mental faculties.
Do not leave the device and its accessories unattended in places accessible to children and/or the disabled.
 8. Do not leave the power supply unattended in places accessible to children and/or people not in full possession of their mental faculties as there is a risk of strangulation with the electric cable.
 9. The medical device may come into contact with the patient via the nebuliser / masks / mouthpiece / nosepiece, components compliant with the requirements of regulation ISO 10993-1, therefore, no allergic reactions and skin irritations may occur.
 10. The product and its parts are biocompatible in accordance with the requirements of regulation EN 60601-1.
 11. Operation of the device is very simple and therefore no further explanations are required other than those indicated in the following user manual.
 12. **ATTENTION:** Do not modify this device without the authorisation of Manufacturer. None of the electrical or mechanical parts of the device have been designed to be repaired by the user. Failure to comply with the above can compromise the safety of the device and, therefore, of its use.
 13. Use of the device in environmental conditions other than those specified in this manual may jeopardise the safety and performance of the device.
 14. The materials in contact with the drugs are thermoplastic polymers with high stability and chemical resistance, widely used in the medical field. However, given the variety and continuous evolution of the used drugs, it is not possible to exclude interactions and, therefore, it is suggested to:
 - Always consume the drug as quickly as possible after its opening
 - Always avoid prolonged contact of the drug with the appropriate tray container (Mesh Chamber) **and immediately carry out the cleaning procedures after each application**
 - Should abnormal situations (e.g. vitrification or cracks) in the medication chamber arise, do not introduce any solution and do not perform inhalation. Contact the technical service specifying the method of use and the type of drug used.
 15. Remember:
 - Do not use the product to nebulise liquids/substances that are not drugs prescribed by your doctor
 - Carry out the aerosol treatment by respecting the proper way to perform it, the doses and combinations indicated by the medical specialist, using only the accessory recommended by them, based on the disease to be treated.

**Do not drop or apply strong shock to the Main Unit or Medicine Cup.
They may be broken and cannot be used or you may suffer electric shock.**

Manufacturer and Distributor cannot be held liable for accidental or indirect damage, if changes have been made to the device, unauthorised repairs and/or technical interventions have been performed or any of its parts have been damaged by accident, improper use and/or abuse.

Any unauthorised intervention, even minimally, on the device voids the warranty immediately and, in any, case does not guarantee correspondence with the technical and safety requirements of the Medical Devices Directive 93/42/EEC (and subsequent changes) and related standards

Do not clean the mesh with cotton balls, brushes or similar! Never touch the mesh. Please follow the sun cleaning instructions indicated in the manual.

PREPARATION AND USE

When first used or always after long time of non-use, we recommend disinfection of the nebulisation chamber, the mask and mouthpiece, so to avoid propagation and infection from miscellaneous bacteria. Then follow the instructions in the Cleaning and Disinfection paragraph of this manual and proceed with the operations described below.

DO NOT INSERT OR NEBULISE substances other than the drugs or medical devices intended for aerosol nebulisation in this device. DO NOT USE WITH Medical Device and/or with dense solutions containing pure hyaluronic acid or dissolved in a high percentage, essential oils, fragrances, etc., or sediment loads (e.g. Pure Thermal Waters) that can block or permanently damage the device. Use with pure hypotonic solutions (e.g. Seawater) can block the device if not properly cleaned immediately afterwards each use.

1 - DRUG INSERTION

FIG A

FIG B

FIG C

A) Open the drug chamber cover by releasing the appropriate safety lever with your fingers.

B) Insert the prescribed drug(s) in compliance with the instructions of your doctor and the drug manufacturer.
ATTENTION: The drug chamber is equipped with a graduated scale indicating the maximum volume of drug that can be inserted for each individual application. This volume is 8ml (8cc). Even if the maximum capacity of the chamber allows up to 12 ml of drug, it is recommended not to exceed the maximum value of 8 ml indicated on the graduated scale.

C) Close the drug chamber cover again, making sure to have properly engaged the safety lever.

Do not drop medication on the Main Unit. If you drop medication, immediately wipe it off with gauze. If you use the unit while it is wet, it may cause trouble or you may suffer electric shock.

On the stop perimeter of the chamber cover, a circular transparent gasket is fixed in its appropriate seat and must never be removed, not even during cleaning operations. Always check that it is in place and in case of accidental loss, do not use the product in order to prevent drug spillage and leakage. If the gasket is lost, you must replace the nebulisation chamber.

2 – PREPARING THE PRODUCT


Once the drug is inserted in the relative chamber, proceed with the following operations:

A) Connect the drug chamber to the main unit making it slide along specific tracks until a click is clearly heard, indicating its proper positioning.

Nebulisation does not start if the drug chamber is not properly engaged. In this case, pressing the ON/OFF button, causes the electronic control of the device to flash, for a few seconds, an BLUE button backlight, that switch off because of missing good connection. If the chamber has been properly engaged, by pressing the ON/OFF button, the button operation BLU light flash slow (cycle of 4 sec.) and permanently until the end of treatment and switch off of device.

B). Connect the accessory indicated by your doctor for the specific therapy to the product. Make sure you use the correct size mask based on the patient (L for adults and S for children), in order to ensure maximum comfort and effectiveness of the treatment.

3 – POWER SUPPLY OF THE DEVICE

KIWI Plus is equipped with electronic power supply management that allows it to function either connected to the rechargeable batteries inside (Original power supply) or by using USB Cable Adapter.

The electronic management of the product automatically identifies the available power source. When the device is connected to the mains electricity supply, it uses this source even if internal batteries are charged, to prevent consuming them.

3.1 Operation with power supply

It is important for the safety of the device and the user, to only use the original power supply that Manufacturer supply as optional accessory. Even the use of power supplies with equal plate data but not of an original type as shown here, is prohibited, because it does not represent an index of equal safety and operational functioning.


Original power supply models:

Shenzhen Huoniu Tecnology Co.,Ltd. - HNEL050100WE, Input:100~240V~, 0.2 AMax, 50/60Hz, Output: 5V $\overline{\text{---}}$, 1A

Shenzhen Huoniu Tecnology Co.,Ltd. - HNB1050100WE, Input:100~240V~, 0.2 AMax, 50/60Hz, Output: 5V $\overline{\text{---}}$, 1A

Shanghai Sucheng Electronic Co.,Ltd. - SC06W-0501000V, Input:100~240V~, 0.15 AMax, 50/60Hz, Output: 5V $\overline{\text{---}}$, 1A

A) The DC connector is located on the bottom of the device. Release the green protective cover as shown in the figure and insert the power supply jack.

B) Connect the power supply to the socket after having ascertained that the type of plug is appropriate. The original power supply provided is of universal voltage type and is, therefore, appropriate for use regardless of the low mains voltage available.

To avoid potentially hazardous conditions, overheating, accidental drops and/or damages, we recommend disconnecting the power supply from the socket immediately after having finished using the device. Do not leave the power supply connected to the socket for a long time.

Do not remove and always close the protective cap (green) of the connector on the bottom when not using the power supply. Always carefully disconnect the power supply from the device to prevent breakage that could affect operation. Any damages resulting from liquid penetration or mechanical failure of the connector are not covered under the Warranty.

DO NOT FORGOT TO SWITCH OFF device at the end of its use. This device DON'T SWITCH OFF AUTOMATICALLY when medicine when medicine is finshed to be nebulised. Only after 20 minutes of continuos operating the device is ALWAYS switch off automatically.

3.2 Operation with internal batteries

The internal lithium battery of this device, rechargeable and not replaceable by the user, allows it to works even when the mains is not available. In full efficiency and 100% charged, the internal battery ensures up to a maximum of 220 minutes of continuous operating of device. With a full charge cycle, the battery is able to power the device to perform up to 28 complete nebulization therapies (1 therapy = 1 total nebulization of 2.5ml of physiological NaCl 0.9%)

The battery is guaranteed for 300 full charge cycles (therefore the useful life of the battery, used and charged in the best conditions, is at least of 1100 hours of operating).

IMPORTANT: The batteries of this device have been pre-charged by the Manufacturer so that device can be used immediately by user. In case instead that device immediately indicates low battery when first switched on (blue light flashing quickly), it is necessary to proceed as soon as possible with a full charge cycle of internal battery.

3.2.1 How to recharge internal battery

A) When the batteries are low (rapid flashing of the BLUE light on the product front side), connect the device to the supplied DC power supply and place the device horizontally on a flat and stable surface, away from heat sources. Then connect the power supply to the power outlet to begin the charging cycle of internal battery.

B) When the power supply is connected to the power outlet, a YELLOW light will light up on the bottom of the device to indicate that the internal battery charging cycle has started regularly. When the charging cycle is complete (battery fully charged) the Yellow light will be replaced by a GREEN light. The charging times changes according to the initial level of charge of the battery. With fully discharged batteries, the full charge time lasts approximately 4 hours (100% charge).

C) The recharging will stop if the device is turned on for nebulization during recharging. Do not stop charging cycle unless urgently needed.

IMPORTANT GENERAL ADVISE:

ALWAYS RECHARGE THE BATTERY WITHIN 1 MONTH FROM THE AUTOMATIC SHUTDOWN FOR LOW CHARGE. IF THE BATTERY IS LEFT COMPLETELY DISCHARGED FOR MORE THAN A MONTH, IT MAY BE DAMAGED PARTIALLY OR TOTALLY


THE BATTERY CAN BE CHARGED FOR A MAXIMUM OF 300 TIMES UNTIL FULLY REACHING ITS FUNCTIONAL SPECIFICATIONS. WHEN MAINTAINED AND USED CORRECTLY, SUBSEQUENT RECHARGES OF BATTERY ARE POSSIBLE BUT IN FUNCTION OF ITS GENERAL INTEGRITY CONDITIONS, THEY MAY PROVIDE LESS AUTONOMY OF USE. ONCE THE BATTERY IS OUT, THE DEVICE CAN ONLY OPERATE WITH THE POWER SUPPLY.

IMPORTANT: For a longer durability and to best exploit battery performance, it is recommended not to carry out partial charges and to recharge the internal battery only when it is discharged (rapid flashing of the front BLUE LED and / or automatic switch-off for low battery).


4 - USING THE DEVICE

Once properly prepared, the device is very easy to use. Just press the ON/OFF button on the front panel to turn the device on or off.

Pressing the button turns on its blue operation backlight and the nebulised medication will start to come out.

At the end of the treatment, when there is no more nebulised medication, press the ON/OFF button to turn the device off.

IMPORTANT: The device can be used continuously, but it is important to remember that it is equipped with a **protective function against accidental start or failed shutdown**, that, **after 20 continuous minutes of operation** (without switch-off), automatically turns the device off.

ATTENTION: During use, the device must be held vertical or tilted towards the patient. The device can operate tilted forwards even up to 90° but if tilted backwards (away from the patient), the drug will be unable to reach the nebulisation Mesh and, therefore, the flow of medication will be interrupted.

To prevent clogging and/or damaging to the mesh, the product must be operated until the drug is fully used up (until steam output has completely stopped). Completely empty the chamber if treatment is stopped before its entire contents have been used. In either case, thoroughly rinse only the drug chamber under warm running water and add a few drops of demineralised water. Engage the chamber to the device and switch the device on for the time required (a few seconds) to drain the last residues. When no more steam is output, switch the product off and clean as described below.

CLEANING AFTER EACH USE

This product must be cleaned after each use, not only for the proper maintenance of the hygienic condition, but to maintain performance and proper operation.

The failed or incorrect cleaning of the drug chamber immediately after its every use, may partially or completely jeopardise its proper operation and consequently require its replacement.

AT THE END OF EACH USE, PROCEED AS FOLLOWS:

Carry out the **SELF-CLEANING MODE** procedure described in the following paragraph, then proceed with the operations described below:



A) After having turned the device off, disconnect the power supply from the socket and from the device itself and store it in a dry and safe place;

B) Pressing the button on the back, release the drug chamber and remove it from the product.

C) Once the main unit has been stored in a safe place (away from sinks or surfaces where it can come into contact with liquids or could fall), open the nebulisation chamber and empty any residual drug;

D) Clean the accessories and drug chamber using water. Carefully rinse the Metal Mesh when using a medication having high surface activity or viscosity such as medication solubilizing agent or expectorant because the medication is often left on the mesh.


Drug chamber, masks and connector, mouthpiece and nosepiece can also be soaked in a 1/1 solution of water and white vinegar for about 1 hour to cleanse these and/or to remove residues and deposits.

DO NOT CLEAN THE MAIN UNIT USING WATER. AVOID CONTACT WITH WET SURFACES AND EXPOSURE TO POSSIBLE SPLASHES
DO NOT CLEAN THE DEVICE, ITS PARTS OR ACCESSORIES IN A DISHWASHER
DO NOT USE DETERGENTS, CHEMICALS OR SCALE REMOVERS TO CLEAN THE DEVICE

E) Place the drug chamber and accessories on a clean and dry cloth and leave to dry completely in the open (about 2 hours) before storing or re-using them.

Do not use cotton balls, brushes or other objects to clean or dry the drug chamber Mesh. **Never touch the mesh!** If necessary, use a dry cotton cloth to dry the accessories and the drug chamber.

F) Use a soft and dry cloth to clean the main unit and the electrodes. Keep the electrodes clean and dry and always remove any residual of medication from the device.

G) Once dry, reassemble the drug chamber to the main door and place them in the box or travel bag, together with the accessories. Keep it all in a clean environment and protected against the weather.

SELF-CLEANING MODE

- Step 1 : Put small amount of clean water in the Medicine Cup.
Step 2 : Turn on the Unit by pressing the Power Button, then hold on the Power Button and don't release for more than 5 seconds, the indicator of Power will be lighted all the time indicating Self-cleaning mode entered.
Step 3: After 30 seconds, the Unit will stop working.
Step 4: Discard the remaining water, and dry the Unit thoroughly.


CARRY OUT THE SELF-CLEANING PROCEDURE AFTER EVERY SINGLE USE

DISINFECTION

Before proceeding with the disinfection and sterilisation, always follow the cleaning procedure described in the previous paragraph.

DRUG CHAMBER, MASKS AND THEIR CONNECTOR, MOUTHPIECE AND NOSEPIECE can be disinfected using denatured alcohol or hypochlorite-based solution, easily available in pharmacies. Before re-using these accessories, rinse them with warm water until all traces of disinfectant have been removed, dry and store in a dry and non-dusty environment

MOUTHPIECE and NOSEPIECE can be disinfected by boiling in water for about **10** minutes or using hot steam sterilising devices (commonly used to disinfect accessories and components for children and infants). Once boiled/disinfecting, leave to cool and dry in the open.

DO NOT USE A MICROWAVE OVEN TO STERILISE THE DEVICE OR ACCESSORIES

MAINTENANCE

The **KIWI Plus** device does not contain any part that requires maintenance and/or lubrication. However, it is necessary to perform a few simple checks to verify the functionality and safety of the device before each use. Pull out the device from the box and **always check** for any visible damage; pay particular attention to cracks in the plastic that could expose some electrical components.

Verify that the drug chamber has not been broken from previous use.

WHEN THE INTERNAL BATTERY WAS AUTO SHUT-OFF BY LOW BATTERY IT HAS TO BE RECHARGED WITHIN ONE MONTH TO AVOID ITS PARTIAL OR FULL DAMAGE

THE INTERNAL BATTERY MUST BE FULLY CHARGED AT LEAST ONCE A YEAR

Upon request, the Manufacturer will provide wiring diagrams, a component list, descriptions, calibration instructions and/or any other information to assist the technical service personnel in repairing the medical device parts.

ACCESSORIES / SPARE PARTS

SEPARATLY AVAILABLE SPARE PARTS

DRUG CHAMBER – REF SP 0279

Soft Touch K+ MASK (L) ADULTS – REF SP 0281

Soft Touch K+ MASK (S) CHILD - REF SP 0282

Soft Touch K+ MASKS CONNECTOR – REF SP 0283

MOUTHPIECE – REF SP 0263

K+ NOSEPIECE – REF SP 0284

AC/DC POWER SUPPLY - REF SP 0262

AC/DC POWER SUPPLY – REF SP 0262/01

EVO Black Carrying Case for KIWI Plus – REF SP 0280

Only use original accessories/spare parts provided and indicated by the

DRUG CHAMBER: For each single patient, we recommend replacing the drug chamber ever 12 months or 350 treatments at the most. In addition, the drug chamber must be replaced in the case of cracks or if the holes in the Vibrating Mesh inside it are partially or totally obstructed by medication, dust, limestone deposits, etc..

In the presence of diseases with risk of infection and microbial contamination, we recommend a personal use of the accessories and drug chamber and compliance with the sterilisation procedures described in this manual (always consult your medical specialist).

Masks and mouthpieces must be replaced if they show visible signs of decay.

Minimum Expected service life: 2 Years for main unit and power supply – MAX 180 operating hours for drug chamber in accordance with the standard test and operation conditions.

Shelf life: max 5 years from the date of manufacture **but max 2 YEARS without application of a full recharging of internal battery**

POSSIBLE PROBLEMS AND SOLUTIONS

Defect	Cause	Remedy
1. Very low or no nebulisation (with device on)	<p>The holes in the vibrating mesh, located inside the drug chamber, are clogged</p> <p>The backlight of Power Button flashed quickly (the interval is less than 0.3 second), and the batteries are low/worn out.</p>	<p>Clean the drug chamber as instructed in the manual. If the problem persists even after cleaning, replace the drug chamber with a new one</p> <p>Recharge the internal battery (full recharging cycle)</p>
2. No nebulisation with device on (power button backlight is on)	<p>The electrodes in the drug chamber or in the main unit are dirty and/or stained</p> <p>Medication or water is collected on the Electrodes of the Main Unit.</p> <p>The Metal Mesh is extremely stained.</p> <p>The Metal Mesh is broken.</p>	<p>Clean the contacts by removing the cause of the poor electric contact.</p> <p>Remove the medication or water.</p> <p>Replace the Medicine Cup with a new one.</p> <p>Replace the Medicine Cup with a new one.</p>
3. Power button backlight is on and after few second turns off but no nebulisation	The drug chamber is not inserted in the main unit or it is not correctly inserted	Mount the drug chamber on the main unit and verify it is correctly locked (limit switch click)
4. Slow nebulisation	Very oily drug	Dilute the drug with a physiological solution
5. No nebulisation and operation backlight of power button is off	The power supply is the wrong type or not correctly inserted in the socket	Verify the use of the original power supply and that it is properly connected to a working wall socket
Defects 1-2-3-4-5	None of the solutions are effective	Contact your dealer or the Manufacturer technical service

The Manufacturer declares that the information contained in this manual is consistent with the technical and safety specifications of the device to which the manual refers. The stated technical specifications are updated to the date of the publication of this document and only apply to the machine to which they are attached. The Manufacturer reserves the right to make changes or improvements without prior notice in this document.

TECHNICAL SPECIFICATIONS

TYPE (93/42/EEC)	Class IIa Medical Device
MODEL	KIWI Plus
POWER SUPPLY	100±240V~ 50/60Hz with provided power supply (mod. HNEL050100WE or SC06W-0501000V) or 3.7V with internal rechargeable Lithium-Iones battery 2600mAh
POWER CONSUMPTION	2.0 W approx.
AUTONOMY WITH INTERNAL BATTERY	approx. 220 minutes (continuously) or 14 days with 2 daily applications(1 application = 1 total nebulization of 2.5ml NaCl0.9%)*
MAX RECHARGING TIME	Approx 4 Hours
USEFUL LIFE OF THE BATTERY	300 cycles of 100% full recharging
RESIDUAL VOLUME (2 spitting)	< 0.1ml
MMAD check production limit 100% of devices (laser diffraction system) ****	< 4.5 µm
MINIMUM NEBULIZATION VOLUME during 100% production checks ****	> 0.25 ml/min (2ml NaCl 0.9%)
MMAD (EN 13544-1) **	1.51 µm
GSD (EN 13544-1) ***	1.75
AEROSOL OUTPUT (EN 13544-1)	1.58 ml
AEROSOL OUTPUT RATE (EN 13544-1)	0.22 (ml/min) 1 min
NEBULIZATION VOLUME	0.30 ml/min (2ml NaCl 0.9%)
INSULATION CLASS	Class II (if used with HNEL050100WE or SC06W-0501000V power supply)
	Internally Powered Equipment (if used with internal batteries)
WEIGHT	Approx.150 g (without accessories)
SIZE	48 (W) x 127 (H) x 58 (D) mm
OPERATION	Continuous
MAX SOUND LEVEL	38.2 dB (A)
MAX MEDICATION CAPACITY	8 ml
OPERATING CONDITIONS	Room temperature: 10 ÷ 40 °C Relative humidity: 30 ÷ 85 %RH (non-condensing) Atmospheric pressure: 700 - 1060 hPa
STORAGE CONDITIONS AND OF TRANSPORT	Room temperature: -20 ÷ 50 °C Relative humidity: 30 ÷ 80 %RH (non-condensing) Atmospheric pressure: 795 - 1013 hPa


* With new battery, intact and with 100% of charge available

** MMAD = Mass Median Aerodynamic Diameter (measured with Marple Personal Cascade Impacter)

*** GSD = Geometric Standard Deviation

**** The parameters shown are used in production during 100% verification of the devices as acceptance limits and therefore represent the minimum threshold value to ensure compliance with the parameters measured in accordance with EN 13544-1

PARTICLE DISTRIBUTION GRAPH AS IN EN 13544-1

N.B.: The measurements and the curves are not valid in case of drugs provided in suspension with high viscosity

MANUFACTURER WARRANTY

CA-MI Srl guarantees its products for a period of at least 24 months from the purchase date. Based on this warranty, CA-MI srl only undertakes to repair or replace the product or its parts that are found to be defective after verification at our premises by the Technical Assistance Service.

The product must be returned accompanied by a description of the defect found and it must be sent together with the Receipt. We always ask for the defect to be specified, in order to perform repair in the shortest possible time. The warranty, with the exclusion of liability for direct or indirect damage, is limited only to material or workmanship defects and is not valid when the pieces returned turn out to have been disassembled, altered or repaired outside the factory or by authorised service centres. The returned device, even if under warranty, must be shipped **CARRIAGE PAID**. Goods always travel at the risk of the buyer, without any liability on the part of CA-MI srl for damage in transit or loss by the carrier, even if sent carriage paid.

The batteries included in this device at the time of purchase are of a commercial type and their charging features may have deteriorated depending on their date of manufacture, their storage and time in the warehouse, therefore, they are not covered by the standard Manufacturer's Warranty.

SYMBOLS USED

	Mark of conformity with Directive 93/42/EEC and subsequent amendments		
	Attention!		
	Consult the user manual		
	Store in a cool dry place		
	Storage temperature		
	Device with BF type applied part (masks/mouthpiece)		
	ON - OFF		
	Device with Class of Insulation II (if used with the provided AC/DC power supply)		
	Alternating current		
	Direct Current		
<p data-bbox="135 732 210 768">IPX2</p>	Degree of protection of an electrical device against accidental or intentional penetration of foreign solid bodies, and protection against the penetration of liquids. <table border="1" data-bbox="267 756 847 852" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%; padding: 5px;"> 1st DIGIT Penetration of solids Not Protected </td> <td style="width: 50%; padding: 5px;"> 2nd DIGIT Penetration of liquids Protected against drops of water diverted up to 15° inclination </td> </tr> </table>	1st DIGIT Penetration of solids Not Protected	2nd DIGIT Penetration of liquids Protected against drops of water diverted up to 15° inclination
1st DIGIT Penetration of solids Not Protected	2nd DIGIT Penetration of liquids Protected against drops of water diverted up to 15° inclination		
	Production lot		
	Serial number		
	Product ID code		
	Manufacturer		

DISPOSAL OF SPENT BATTERIES (Directive 2006/66/EC)

This symbol on the product indicates that the batteries should not be treated as household waste. Ensuring that the batteries are correctly disposed of contributes to preventing potential negative consequences for health and the environment, which could otherwise be caused by their inappropriate disposal. Recycling materials helps to preserve our natural resources. Take the spent batteries to collection points for recycling. For more information about collection and recycling of spent batteries or the product, please contact your Municipality, your local waste disposal service or the shop where you purchased the device.

WARNING REGARDING THE CORRECT DISPOSAL OF THE PRODUCT ACCORDING TO THE EUROPEAN DIRECTIVE 2012/19/EU-WEEE:

After use, the product must not be disposed of with household waste. It can be delivered to specific separate collection centres provided by the municipal administration, or returned to the distributor when purchasing a new device of the same type and having the same functions. Disposing of the product separately avoids possible negative effects on the environment and human health and allows the recovery of materials of which it is composed in order to obtain significant savings in energy and resources. The symbol on the data label indicates the separate collection of electrical and electronic equipment.

Attention: An incorrect disposal of electrical and electronic equipment could lead to fines.

RISK OF ELECTROMAGNETIC INTERFERENCE AND POSSIBLE REMEDIES

This section contains information on device conformity with EN 60601-1-2. The KIWI Plus aerosoltherapy device is an electromedical device that requires special precautions regarding electromagnetic compatibility and needs to be installed and commissioned according to the information specified in the accompanying documents. Portable and mobile radio communication devices (mobile phones, transceivers, etc.) may interfere with the medical device and should not be used in close proximity with, adjacent to or on top of the medical device. If such use is necessary and unavoidable, special precautions should be taken so that the electro-medical device functions properly in its intended operating configuration (for example, constantly and visually checking for the absence of anomalies or malfunctions). The following tables supply information regarding the EMC (Electromagnetic Compatibility) characteristics of the electro-medical device.

Manufacturer's guide and declaration - Electromagnetic Emission		
The KIWI Plus aerosol is for use in the below specified electromagnetic environment. The Customer and/or user of the KIWI Plus aerosol must ensure that the device is used in such an environment		
Emissions Test	Conformity	Guide to the electromagnetic environment
CISPR11 Radiated / Conducted Emissions	Unit 1	The KIWI Plus aerosol uses RF energy only for its internal operation. Therefore, its RF emissions are very low and do not cause any interference with any nearby electronic device.
CISPR11 Radiated / Conducted Emissions	Class [B]	The KIWI Plus aerosol is suitable for use in all environments, including domestic, and those directly connected to the public distribution network providing power to environments used for domestic purposes.
Harmonic EN 61000-3-2	Class [A]	
Voltage fluctuations / flicker EN 61000-3-3	Compliant	

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