iX Auto CPAP System

Continuous Positive Airway Pressure User Manual

Be well, Live well
Wellell

CONTENTS

1. USER'S MANUAL	1
2. MANUAL DE INSTRUCCIONES	45
3. INSTRUÇÕES DE UTILIZAÇÃO	92
4. INSTRUCTIONS D'UTILISATION	139
5. GEBRAUCHSANLEITUNG	188
6. GEBRUIKERSHANDLEIDING	237
7. MANUALE DI ISTRUZIONI	282

MODEL NO.: 9S-010600 & 9S-010620
PLEASE READ THE MANUAL BEFORE USE

IMPORTANT SAFEGUARDS – PLEASE MAKE THIS PRECAUTIONS

LIST AVAILABLE TO ALL USERS AND PROVIDERS

READ ALL INSTRUCTIONS BEFORE USING THE DEVICE

NOTE, CAUTION, AND WARNING STATEMENTS

-  **NOTE:** Indicates information that you should pay special attention to.
-  **CAUTION:** Indicates correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.
-  **WARNING:** Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.

DANGER – To reduce the risk of electrocution

- (1). **DO NOT** place or store product where it can fall into or be pulled into a tub or sink.
- (2). **DO NOT** place or immerse the unit into water or other liquids.

WARNINGS – To reduce the risk of burns, electrocution, fire or injury to the user or operator.

- (1). THIS DEVICE IS NOT INTENDED FOR LIFE SUPPORT.
- (2). The airflow for breathing generated by this device may be as much as 5°C (41°F) higher than the room temperature. This device should not be used if the room temperature is warmer than 35°C (95°F) to prevent the airflow temperature from exceeding 43° C (109° F) and causing irritation to your airway.
- (3). This machine should be used only with masks (and connectors) recommended by the manufacturer, or by your physician or respiratory therapist.
- (4). At low CPAP pressure, some exhaled gas may remain in the mask and be re-breathed. Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocation.
- (5). If you are using a full face mask (a mask covering both your mouth and your nose), the mask must be equipped with a safety (Anti-asphyxia) valve.
- (6). Do not place the heater base and the humidifying chamber higher than the patient's position, otherwise, condensed water droplets may enter into the patient's airway.

- (7). Equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with nitrous oxide or oxygen.
- (8). Close supervision is necessary when this product is used by, on, or near children or invalids. Choking may result from a child swallowing a small part that has become detached from the device or its accessories.
- (9). Use this product only for its intended purpose as described in this manual. DO NOT use auxiliary attachments that are not approved by the manufacturer.
- (10). Use only Wellell or confirmed by Wellell parts and accessories with the device. Non-Wellell parts may reduce the effectiveness of the treatment, potential allergic reactions or damage the device.
- (11). Never operate this product if it has a damaged cord or plug. If it is not working properly, or has been dropped, damaged, or immersed in water, return the product to a service center for examination and repair.
- (12). Never block the air outlets on this product. DO NOT place the unit on a soft surface such as a bed or couch, where the air outlets may become obstructed. Keep the air outlets free from lint, hair, and other obstructive particles.
- (13). Never drop or insert any object into any opening on the unit or into the tubing.
- (14). When disposing the unit because of damage or failure, you must comply with the regulations of your particular region or country for disposing of such equipment.
- (15). No modification of this equipment is allowed. If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.
- (16). Do not leave long lengths of air tubing around the top of your bed. It could twist around your head or neck while you are sleeping and lead to strangulation.
- (17). Make sure the environment around the machine is dry and clean. Dust and foreign particles may affect the treatment. Keep the air inlet on the back of the machine clear to prevent overheating and damage of the device. Do not place the machine near a source of hot or cold air. Extreme cold or hot environment may damage the user's respiratory airway.
- (18). If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance between devices or turn off the mobile phone.
- (19). Do not connect the device to the personal computer for data downloading during the treatment. This may cause the CPAP device failure.
- (20). U.S. Federal law restricts this device to sale by or on the order of a licensed physician.
- (21). Make sure the proper placement and positioning of the MASK on the face is critical to the consistent operation of this equipment.

- (22). Medical electrical equipment needs special precautions regarding EMC and needs to be installed according to EMC information. Contact your home care provider regarding EMC installation information.
- (23). Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- (24). Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- (25). Sources of oxygen must be located more than 1 m from the equipment to avoid the risk of fire and burns.
- (26). Failure to use a mask or accessory that minimizes rebreathing of carbon dioxide or permits spontaneous breathing can cause asphyxiation.

CAUTION – Reducing the risk of damage to the device.

- (1). Mobile RF communications equipment can affect medical electrical equipment.
- (2). Ensure that the SD card sealed cover and filter back cover are both closed whenever any of the accessories are not installed.
- (3). Do not use any extension cords with this device.
- (4). Make sure the filter area at the back of the device is not blocked and check and change the inlet filter for each specific time (The detail as Chapter 6 Cleaning & Maintenance).
- (5). Do not install a wet filter into the device.
- (6). Do not place the device directly onto carpet, fabric, or other flammable materials.
- (7). Do not place the device in or on any container that can collect or hold water.
- (8). Tobacco smoke may cause tar build-up within the device, which may result in the device malfunctioning.
- (9). Make sure that the DC power cord securely insert into your therapy device before using.
- (10). Only use Wellell approved DC Power Cord and Adapter. Use of any other system may cause damage to the device.
- (11). The equipment must not be covered or positioned in such a way that the operation or performance of the equipment is adversely affected.
- (12). Be careful not to place the device where it can be taken or touched by children and pets.
- (13). Place the device away from sources of heat and moisture to avoid damaging the equipment.

1. Introduction

This manual should be used for initial set up of the system and saved for reference purpose.


1.1. General Information

Obstructive Sleep Apnea (OSA) is a condition that an intermitted and repetitive obstruction of the upper respiratory tract causes a complete (apnea) or partial (hypopnea) block of breathing airflow during sleep. The syndrome varies depending on the degree of relaxation of the tongue and soft palate muscle.

The most common treatment for OSA is Continuous Positive Airway Pressure (CPAP). CPAP devices can deliver a constant air pressure into your upper airway via a mask. This constant air pressure can keep your airway open during sleep, and therefore prevents the OSA.

This device is a micro-processor controlled continuous positive airway pressure device. It features the illuminated, menu-driven LCD display, universal power supply and ramp time adjustment. The ramp time adjustment and ultra-quiet operation ensure you to fall asleep comfortably while air pressure slowly builds up to treatment level. The user compliance meter records the total system's operating time for physician's reference.

The system has been tested and successfully approved to the following standards:

 **For US and CANADA only**

E348970
53DG

Medical Equipment- CPAP

WITH RESPECT TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH

ANSI / AAMI ES60601-1 (2005) + AMD 1 (2012)

CAN / CSA-C22.2 No. 60601-1:14

 **The below description for is for Canada only: The product has been tested with medical equipment and meets ANSI / AAMI ES60601-1 (2005) + AMD 1 (2012) standards; CAN / CSA-C22.2 No. 60601-1: 14. Preventing electric shock, fire and the risk of physical injury.**

1.2. EMC Warning Statement

This equipment has been tested and found to comply with the limits for medical devices to the

IEC/EN 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Connect the equipment to an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service technician for help.

1.3. Intended Use

This device is intended to provide continuous positive airway pressure (CPAP) for the treatment of adult Obstructive Sleep Apnea (OSA) in spontaneously breathing patients weighing over 30 kg (66 lbs). It is for single patient reuse in the home environment.

Caution: Some patients might have pre-existing contraindications for CPAP therapy, or might experience some potential side effects of using a CPAP device. Please consult your doctor if you have any questions concerning your therapy.

Warning: Equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

1.4. Contraindications

Studies have shown that the following pre-existing conditions may contraindicate the use of CPAP therapy for some patients. Should you have any of these conditions, your physician will determine if CPAP therapy is appropriate for you.

- (1). Uncooperative or extremely anxious patient
- (2). Reduced consciousness and inability to protect their airway
- (3). Unstable cardiorespiratory status or respiratory arrest
- (4). Trauma or burns involving the face
- (5). Facial, esophageal, or gastric surgery
- (6). Air leak syndrome (pneumothorax or bronchopleural fistula)
- (7). Copious respiratory secretions
- (8). Severe nausea with vomiting
- (9). Severe air trapping diseases with hypercarbia asthma or chronic obstructive pulmonary disease (COPD)
- (10). Pathologically Low Blood Pressure

1.5. Adverse effects

The following side effects may be noticed during the CPAP therapy with the device:

- (1). Dry or stuffy nose
- (2). Dry mouth or throat
- (3). Nosebleed
- (4). Bloating
- (5). Ear or sinus discomfort
- (6). Eye irritation
- (7). Skin rashes

If unusual chest pain, severe headache, or increased breathlessness happens, patients should report to their prescribing physician. An acute upper respiratory tract infection may require temporary discontinuation of the treatment.

2. Product Description

2.1. System Contents

Your iX Auto system may include the following items:

- (1). Device
- (2). Adapter
- (3). AC power cord
- (4). Light grey tubing –15mm
- (5). Carrying bag
- (6). Ultra-Filter
- (7). User manual
- (8). SD card
- (9). Micro USB cable (for Technician)

Optional Items:

- (1). Heated Humidifier with Water Chamber
- (2). Heated Tubing

NOTE: Please use a SD card which complies with SDHC standard. Before using it, please format it to FAT 32 to ensure correct data collection.

NOTE: The content of standard package might vary or change without advanced notification

CAUTION: Patient should not connect the device to a personal computer for data downloading. This may cause the CPAP device failure.

2.2. System Overview

- (1). LCD Screen
- (2). Power inlet
- (3). Filter cover
- (4). Air outlet
- (5). SD card cover
- (6). Humidifier (optional)
- (7). Water chamber (optional)
- (8). Air outlet
- (9). SD card slot
- (10). Micro-USB cable slot (for Technician only)
- (11). Release button of humidifier

3. Installation

3.1. Unpacking

To secure its contents inside, the iX Auto CPAP device and accessories are bundled in a paper-packaged box. Unpack this box by removing the device and its accessory and checking for any damage, which may have occurred during shipping. If there are damages, please contact your care provider immediately.

3.2. Setting Up

3.2.1. Using device only

Place the device on a stable level surface. Plug the power connector into the rear of the device. Connect one end of the power cord into the power supply unit and the other end into the power outlet.

Connect the Light grey tubing (or optional heated tubing) firmly to the air outlet located on the side of the device.

Caution: Connect the sensor port of heated tubing with device.


Connect the free end of the air tubing firmly onto the assembled mask. See the mask user guide for detailed information.

Recommended masks are available on www.wellell.com

3.2.2. Using with humidifier and heated tubing (Optional)

iX Auto CPAP device can be used with iX Heated Humidifier (Optional) WITH Light grey tubing or heated tubing.

	<p>Place the heated humidifier on a stable level surface. Release the bottom pad from above without pressing the button. The heated humidifier system may reduce nasal dryness and irritation by providing adequate moisture and heat to the airflow.</p>
	<p>Connect the device with heated humidifier horizontally.</p> <p> NOTE: Please make sure the humidifier is inserted to device completely to reach desired performance.</p>
	<p>Hold the water chamber at the top and bottom and take the chamber out by the right side.</p> <p> Caution: Beware of burns when removing the water chamber.</p>
	<p>Open the top cover and fill it with distilled water up to the maximum water level mark, then close it. Slightly wipe the exterior of the water chamber with a dry cloth. Insert the chamber back in the humidifier and make sure both are properly connected.</p> <p> Warning: Do not overfill the water chamber as water may enter the device.</p> <p> Caution: After filling the water, keep the water chamber at level position before assembling with the device.</p> <p> Caution: Please fill the water chamber only with room temperature pure water (about 20°C/68°F).</p>

	<p>After assembling th humidifier with the device, the user can connect the Light grey tubing (or heated tubing) with the rotatable air outlet located on the top of the humidifier.</p> <p> Caution: Connect the senser port of heated tubing with device.</p> <p> Warning: Forcing the connection might cause damage to the connector either on the humidifier or tubing.</p>
	<p>Plug the power connector into the rear of the device. Connect one end of the power cord into the power supply unit and the other end into the power outlet.</p>
	<p>Connect the free end of the air tubing firmly onto the assembled mask. See the mask user guide for detailed information.</p> <p>Recommended masks are available on www.wellell.com</p>

3.2.3. Disconnect accessories – Humidifier, Light grey tubing and Heated tubing

	<p>Humidifier</p> <p>Slightly press the release button on the humidifier, and the device will disconnect from humidifier. Now you can horizontally pull the device away from the humidifier.</p> <p> Warning: Always empty the water inside the water chamber of humidifier after every usage.</p> <p> NOTE: To avoid damage, please do not forcefully disconnect the device without pressing the release button.</p>
	<p>Light grey tubing</p> <p>Tubing could be horizontally removed from the device by holding both sides of the cuff.</p> <p> NOTE: To avoid damage, please do not hold the conjunction area between the Light grey tubing and cuff when removing the tubing.</p>
	<p>Heated Tubing (device only)</p> <p>Heated tubing could be horizontally removed from device by pressing both sides of the cuff.</p> <p> NOTE: To avoid damage, please do not hold the conjunction area between the Tubing and cuff when removing the Tubing.</p>
	<p>Heated Tubing (with humidifier)</p> <p>Heated Tubing could be horizontally removed from humidifier by pressing both sides of the cuff.</p> <p> NOTE: To avoid damage, please do not hold the conjunction area between the Tubing and cuff when removing the Tubing</p>

4. Operation

 **NOTE:** Always read the operating instructions before use.

4.1. Control Panel Description

	<p>1. START / STANDBY Button Start treatment : press "START / STANDBY" (screen : show therapy pressure "XX.X cmH₂O") Stop treatment : press "START / STANDBY" (screen : show current time "XX:XX")</p>
	<p>2. MENU Button Enter main menu: (in standby mode) press "MENU" to enter main menu. There are three settings in the main menu: Comfort, Setting and Report. Return main menu: (in sub-menu) press "MENU" to return to main menu. (or you can press "Back" to return to main menu)</p>
	<p>3. DIAL Select / Conform button (for main menu): rotate "Dial" to select desire menu, then press "Dial" to enter sub-menu. (or you can select "Back" for returning to Standby Mode) Select / Conform button (for sub-menu): rotate "Dial" to select</p>

	<p>desire sub-menu, then press "Dial" to start adjusting parameters. (Please refer to 4.5 for more details)</p> <p style="text-align: center;"> NOTE: In therapy mode, the "Dial" could be used as Quick Start button for adjusting the level of heated Tubing and humidifier if the device has connected with an optional accessory. (Please refer to 5)</p> <p>Level adjustment of humidifier: rotate "Dial" to adjust level, then press "Dial" to confirm it.</p> <p>Level adjustment of heated Tubing: rotate "Dial" to adjust level, then press "Dial" to confirm it.</p>
--	---

4.2. Getting Started – First time use

For first-time use, follow the steps below:

- (1). Make sure that your doctor or care provider has initialized and configured the appropriate treatment values and the correct date and time into the unit. Please make a consult with your doctor periodically to reassess the setting of the therapy for effectiveness.
- (2). Make sure that you have a mask with headgear. If not, your doctor or care provider can recommend the type of mask and headgear you must obtain.
- (3). Attach power cord to an electrical outlet (please refer to 3.2). The device will turn on automatically
- (4). Attach tubing and mask (please refer to 3.2).
- (5). Use the **"MENU"** button and **"Dial"** to set up Comfort, Setting and to read Report.

 **NOTE: these values may have already been set by your doctor or care provider.**

 **NOTE: make sure SD has been inserted into the SD slot at the back side of the main device.**


- (6). You can now place the Headgear and Mask over your head and face.
- (7). Press the **"START / STANDBY"** button – Airflow to the mask will begin. (if the **"Auto ON/OFF"** has been turned on, airflow will automatically begin when you place the mask on your face properly.)

 **NOTE: If power is interrupted during therapy, the device will return to the Home screen once power is restored. You may resume therapy as needed.**

4.3. Function Description

Stand-by




Function menu

⚠ **NOTE:** Rotate and press the "Dial" to enter the icon , the screen will return to stand-by mode.


4.4. Comfort

Screen	Descriptions
	<p>Ramp Time</p> <p>The Ramp Time feature allows the patient to fall sleep at a lower, more comfortable pressure, after which the unit will gradually ramp up to the appropriate therapy pressure. This function is the time (in minutes) it will take for the unit to “ramp” up to the appropriate therapy pressure.</p> <p>Setting: Rotate "Dial" to highlight Ramp Time → press "Dial" → rotate "Dial" to adjust time → press "Dial"</p> <p>(Options: OFF, 5, 10, 15, 20, 25, 30, 35, 40 and 45 minutes; increment: 5mins)</p> <hr/> <p>Auto ON</p> <p>The device will start treatment automatically once the patient places the headgear and mask over the head and face when turning on this function.</p> <p>Setting: Rotate "Dial" to highlight Auto ON → press "Dial" → rotate "Dial" to set ON or OFF this function. → press "Dial"</p> <p>⚠ NOTE: Auto ON function can only be activated in stand-by mode.</p>

Screen	Descriptions
	<p>Auto OFF</p> <p>The device will stop treatment automatically once the patient removes the mask.</p> <p>Setting: Rotate "Dial" to highlight Auto OFF → press "Dial" → rotate "Dial" to set ON or OFF this function. → press "Dial"</p> <hr/> <p>Humidity Level (optional)</p> <p>Setting (in Standby mode): Rotate "Dial" to highlight Humidity Level → press "Dial" → rotate "Dial" to adjust level → press "Dial"</p> <p>Setting (in Therapy mode): please refer to 5 Quick Start Guide for Heated Humidifier and Heated Tubing (Optional).</p> <p>(Options: OFF, 1~6 levels)</p> <p>Level 1: provides lower humidity</p> <p>Level 6: provides higher humidity</p> <p style="text-align: right;"> NOTE: Only when Heated Humidifier has been attached.</p>
	<p>Heated Tubing Setting (optional)</p> <p>Setting (in Standby mode): Rotate "Dial" to highlight Heated Tubing Setting → press "Dial" → rotate "Dial" to adjust level → press "Dial"</p> <p>Setting (in Therapy mode): please refer to 5 Quick Start Guide for Heated Humidifier and Heated Tubing (Optional)</p> <p>(Options: OFF, 1~6 levels)</p> <p>Level 1: provides lower temperature</p> <p>Level 6: provides higher temperature</p> <p style="text-align: right;"> NOTE: Only when Heated Tubing has been attached.</p> <hr/> <p>Mask Type</p> <p>Setting: Rotate "Dial" to highlight "Mask Type" → press "Dial" → rotate "Dial" to select mask type →</p>

Screen	Descriptions
	<p>press "Dial"</p> <p>(option: Nasal, Pillow, or Full face)</p> <p> NOTE: To enhance the treatment performance, please make sure make type and mask type settings are the same.</p>
	<p>Run Mask Fit</p> <p>Setting: Rotate "Dial" to highlight Run Mask Fit → press "Dial" → (entering testing screen)</p> <p>Screen:</p> <p>Good: test pass</p> <p>Please try again: test fail</p> <p> NOTE: If Auto ON set "On", please disconnect the mask and tubing before able this function. When press Dial to start "run mask it", and then quickly connect the mask and tubing.</p> <p> NOTE: make sure that mask type is correctly selected before running mask fit.</p>
	<p>Back</p> <p>Setting: Rotate "Dial" to highlight Back → press "Dial" (to return to main menu screen)</p>

4.5. Settings

Screen	Descriptions								
<div style="background-color: #003366; color: white; padding: 5px; text-align: center; font-weight: bold;">Setting</div> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="padding: 2px;">Airplane Mode</td> <td style="text-align: right; padding: 2px;">Off</td> </tr> <tr> <td style="padding: 2px;">Wi-Fi</td> <td style="text-align: right; padding: 2px;">On</td> </tr> <tr> <td style="padding: 2px;">Cellular</td> <td style="text-align: right; padding: 2px;">On</td> </tr> <tr> <td style="padding: 2px;">Backlight</td> <td></td> </tr> </table>	Airplane Mode	Off	Wi-Fi	On	Cellular	On	Backlight		<p>Airplane mode</p> <p>Setting: Rotate "Dial" to highlight Airplane Mode → press "Dial" → rotate "Dial" to turn ON or turn OFF this function. (Wi-Fi will be turned off at the same time) → press "Dial"</p> <p>Screen: will be shown on the screen when turning on this function. </p>
	Airplane Mode	Off							
	Wi-Fi	On							
Cellular	On								
Backlight									
<p>Wi-Fi</p> <p>Setting: Rotate "Dial" to highlight Wi-Fi → press "Dial" → rotate "Dial" to turn ON or turn OFF this function. (please make sure Airplane Mode is turned off) → press "Dial"</p> <p>Screen:</p> <p> Wi-Fi is connected.</p> <p> Wi-Fi is not connected (or the Wi-Fi signal is too weak).</p>									
<p>Cellular (optional)</p> <p>Setting: Rotate "Dial" to highlight Cellular → press "Dial" → rotate "Dial" to turn ON or turn OFF this function. (please make sure Airplane Mode is turned off) → press "Dial"</p> <p>Screen:</p> <p> Cellular is connected.</p> <p> Cellular is not connected (or Cellular is too weak).</p> <p>Data will be transmitted automatically when each therapy finishes. Do not unplug power cord during data transmission. </p> <p style="text-align: right;"> NOTE: Only for iX series 9S-010620</p>									

Screen	Descriptions								
	<p>Backlight</p> <p>Setting: Rotate "Dial" to highlight Backlight → press "Dial" (for entering next page)→ rotate "Dial" to adjust backlight → press "Dial"</p>								
<div style="background-color: #000080; color: white; padding: 5px; text-align: center; font-weight: bold;">Setting</div> <table border="0" style="width: 100%; background-color: #333; color: white;"> <tr> <td style="padding: 5px;">Alarm Clock Setting</td> <td style="padding: 5px; text-align: right;">On</td> </tr> <tr> <td style="padding: 5px;">Alarm Clock Time</td> <td style="padding: 5px; text-align: right;">13:30</td> </tr> <tr> <td style="padding: 5px;">Reminders</td> <td></td> </tr> <tr> <td style="padding: 5px;">Language</td> <td style="padding: 5px; text-align: right;">English</td> </tr> </table>	Alarm Clock Setting	On	Alarm Clock Time	13:30	Reminders		Language	English	<p>Alarm Clock Setting</p> <p>Setting: Rotate "Dial" to highlight alarm clock → press "Dial" → rotate "Dial" to turn ON or turn OFF this function.</p> <p>Screen:  will be shown on the screen when turning on this function. Press any buttons to mute the audio alarm.</p> <p>Alarm Clock Time</p> <p>Setting: Rotate "Dial" to highlight Alarm Clock Time → press "Dial" → rotate "Dial" to set "Hours" → rotate "Dial" to set "Minutes" → press "Dial"</p> <p>Reminders</p> <p>Setting: Rotate "Dial" to highlight Reminders→ press "Dial" to view the reminder of Filter, Water Chamber, Tube replacement.</p> <p>Languages</p> <p>Setting: Rotate "Dial" to highlight Language → press "Dial" → rotate "Dial" to set the preferred language → press "Dial"</p> <p>(Options: English, Español, 日本語, Nederlands, Deutsch, Français, Italiano, Português, 繁體中文, 简体中文)</p>
Alarm Clock Setting	On								
Alarm Clock Time	13:30								
Reminders									
Language	English								
<div style="background-color: #000080; color: white; padding: 5px; text-align: center; font-weight: bold;">Setting</div> <table border="0" style="width: 100%; background-color: #333; color: white;"> <tr> <td style="padding: 5px;">About</td> </tr> <tr> <td style="padding: 5px;">Back</td> </tr> </table>	About	Back	<p>About</p> <p>Screen: The updated date and software version are displayed.</p> <p>Back</p> <p>Setting: Rotate "Dial" to highlight Back → press "Dial" (to return to main menu screen)</p>						
About									
Back									

4.6. Report (Compliance information)

Rotate **"Dial"** to set the period of time. You can see the compliance data in the setting period. The compliance data can only view. Your home care provider may periodically ask you for this information.


Screen	Descriptions								
<div style="background-color: #003366; color: white; padding: 5px; text-align: center; font-weight: bold;">Report</div> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="padding: 5px;">Period</td> <td style="text-align: right; padding: 5px;">1 Day</td> </tr> <tr> <td style="padding: 5px;">Day Used</td> <td style="text-align: right; padding: 5px;">0 Days</td> </tr> <tr> <td style="padding: 5px;">Usage Hour</td> <td style="text-align: right; padding: 5px;">0.0 Hrs</td> </tr> <tr> <td style="padding: 5px;">Avg. Usage Hour</td> <td style="text-align: right; padding: 5px;">0.0 Hrs</td> </tr> </table>	Period	1 Day	Day Used	0 Days	Usage Hour	0.0 Hrs	Avg. Usage Hour	0.0 Hrs	<p>Period Set the time interval covered by the <i>"Report"</i> (options: 1 Day, 7 Days, 30 Days, 90 Days)</p> <hr/> <p>Day Used Total treatment days during the selected period.</p> <hr/> <p>Usage Hour Total treatment hours during the selected period.</p> <hr/> <p>Avg. Usage Hour Average treatment hours per day during the used day.</p>
Period	1 Day								
Day Used	0 Days								
Usage Hour	0.0 Hrs								
Avg. Usage Hour	0.0 Hrs								
<div style="background-color: #003366; color: white; padding: 5px; text-align: center; font-weight: bold;">Report</div> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="padding: 5px;">Compliance Rate</td> <td style="text-align: right; padding: 5px;">0%</td> </tr> <tr> <td style="padding: 5px;">AHI</td> <td style="text-align: right; padding: 5px;">0.0</td> </tr> <tr> <td style="padding: 5px;">Avg. Pressure</td> <td style="text-align: right; padding: 5px;">0.0 cmH₂O</td> </tr> <tr> <td style="padding: 5px;">90%</td> <td style="text-align: right; padding: 5px;">0.0 cmH₂O</td> </tr> </table>	Compliance Rate	0%	AHI	0.0	Avg. Pressure	0.0 cmH ₂ O	90%	0.0 cmH ₂ O	<p>Compliance Rate The percentage of the day of treatment time more than 4 hours per day during the selected period.</p> <hr/> <p>AHI The number of apnea plus the number of hypopneas on average during per hour of sleep.</p> <hr/> <p>Average Pressure Average pressure during the selected period.</p> <hr/> <p>90% The average pressure at or below 90% of the treatment period.</p>
Compliance Rate	0%								
AHI	0.0								
Avg. Pressure	0.0 cmH ₂ O								
90%	0.0 cmH ₂ O								
<div style="background-color: #003366; color: white; padding: 5px; text-align: center; font-weight: bold;">Report</div> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="padding: 5px;">Leakage</td> <td style="text-align: right; padding: 5px;">0.0 L/min</td> </tr> <tr> <td style="padding: 5px;">Check Code</td> <td style="text-align: right; padding: 5px;">000</td> </tr> <tr> <td style="padding: 5px;">Back</td> <td></td> </tr> </table>	Leakage	0.0 L/min	Check Code	000	Back		<p>Leakage The average leakage is detected during the treatment period.</p> <hr/> <p>Check code It the replacement of phone-in data. Your care provider might ask you about this data regularly.</p> <hr/> <p>Back Setting: Rotate "Dial" to highlight Back → press "Dial" (to return to main menu screen)</p>		
Leakage	0.0 L/min								
Check Code	000								
Back									

4.7. Turn Off the Device

- (1). Press the **“START / STANDBY”** button – Device will be switched to Standby mode.
- (2). Disconnect and remove the power cord from the device.
- (3). Disconnect all accessories such as Light grey tubing, heated Tubing and humidifier (please refer to 3.2.3 for more details)
- (4). Clean the device, tubing and accessory regularly (please refer to 6 for more details).

5. Quick Start Guide for Heated Humidifier and Heated Tubing (Optional)

To enhance operating efficiency and convenience, iX Auto CPAP device build in two quick start functions during therapy mode.

	<p>Below two quick start icon will be shown on the screen when you attach the heated Tubing and humidifier.</p>  <p>Setting (in Therapy mode): Rotate “Dial” to highlight Heated Tubing Setting or Humidifier icon → press “Dial” → rotate “Dial” to adjust level → press “Dial”</p>
---	--

6. Cleaning & Maintenance

6.1. Cleaning the Device

The device should be checked and dusted regularly (at least every 30 days).

- (1). Unplug the device.
- (2). Wipe the outside of device with a cloth dampened with water and mild detergent (please disconnect accessories if necessary. Please refer to 3.2.3).
- (3). Device should be air-dried thoroughly before use. Inspect the device and all circuit parts for any damage after cleaning and replace if necessary.

 **WARNING: To avoid electrical shock, always unplug the power cord from the wall outlet before cleaning the device.**

 **WARNING: Don't try to open this device. Repairs and internal servicing should only be performed by an authorized service agent. Do not drop any object into the air Tubing or air outlet.**

 **WARNING: Please do not repair the device during operation.**

6.2. Cleaning the Humidifier (and water chamber)

The humidifier should be checked and cleaned daily (or after every time use) to prevent the growth of germs that might affect the user's health.

- (1). Disconnect accessories if necessary. (Please refer to 3.2.3)
- (2). Wipe the outside of the humidifier with a cloth dampened with water and mild detergent
- (3). Open the water chamber and discard any remaining water.
- (4). Wash the water chamber with warm water and mild detergent.
- (5). Rinse the water chamber and allow to dry out of direct sunlight and/or heat. Inspect all parts for any damage after cleaning and replace if necessary.

6.3. Cleaning the Tubing

The air tubing should be checked and cleaned daily (or after every time use) to prevent the growth of germs that might affect the user's health (please also refer to the air tubing instruction manual).

- (1). Disconnect the air tubing from the air outlet of the device (Please refer to 3.2.3).
- (2). Wash the air tubing in warm water using mild detergent. Rinse thoroughly, hang, and allow to dry. All items of the-air tubing are subject to normal wear and tear and may eventually be replaced. Replace the air tubing if they are damaged.

 **WARNING: Do not use harsh soaps, detergents, solvents or alcohols to clean the air tubing. These solutions may reduce the life of the product.**

 **WARNING: Do not wash or dry the air tubing at a temperature above 70°C (158°F)**

 **WARNING: Do not use any cleaner containing fragrance or conditioners as they will leave a residue.**

6.4. Air Filter

For an optimal operation of the device, the Ultra filter is disposable and should be replaced after 6 months of use or sooner if it appears dirty; however, depending on the air quality, the replacement time may vary.

- (1). Open back cover and remove the dirty filter from the enclosure on the rear of the device.
- (2). Replace a new filter and close back cover.

 **WARNING: The Ultra filter is not washable or reusable.**

 **WARNING: Only WELLELL-supplied filters should be used as replacement filters.**

 **CAUTION: Dirty air filter may cause high operating temperatures that affect device performance.**

7. Therapy data

iX Auto records your therapy data for you and your care provider so they can view and make changes to your therapy if required. The data is recorded and then transferred to your care provider wirelessly, if a wireless network is available, or via an SD card.

7.1. SD card

Your therapy data can be transferred to your health care provider via the SD card. Your care provider may ask you to send the SD card by mail or to bring it in. When instructed by your care provider, remove the SD card.

To remove the SD card:

- (1). Open the SD card cover.
- (2). Push in the SD card to release it. Remove the SD card from the device. Place the SD card in the protective folder and send it back to your care provider.

 **NOTE: The SD card should not be used for any other purpose.**

 **NOTE: Do not remove the SD card from the device when data is being written to the card.**

7.2. Wireless Technology

Your iX Auto is equipped with a wireless module, so you can transfer your therapy data to SleepWell Track. SleepWell Track is a web-based system designed to help patients enhance their sleep therapy.

How to enable wireless communication:

7.2.1. Wi-Fi

Wi-Fi communication allows your therapy data to be wirelessly transmitted to your health care provider to enhance the quality of your treatment.

- (1). Make sure the Wi-Fi function is on in "Setting" mode. (please refer 4.4) The Wireless signal icon will display at the top of your screen.
- (2). Your device can be used with WELLELL's APP (SleepWell). Please download and install the app from the App Store or Google Play.
- (3). Enable Bluetooth on your mobile device.
- (4). Click on the APP and create account.
- (5). On the Navigator Bar, click "Wifi Setting". Setup the connection to Wi-Fi.

NOTE: The APP will be automatically paired with IX devices via Bluetooth without any configuration.

- (6). Tap a listed Wi-Fi Hotspots
- (7). Enter the correct Wi-Fi network password, and press Connect. This will complete the device's connection to the wireless network.
- (8). The data is automatically transmitted to SleepWell Track after therapy has completed. In order to make sure that your data is transferred, leave your device connected to the power outlet at all times and make sure that the Wi-Fi is on.

NOTE: Please be aware that the availability and quality of the wireless network may be affected by location, buildings, number of users, terminal equipment and other relevant factors.

NOTE: You may turn on/off Wi-Fi function in Wi-Fi setting menu or turn on airplane mode directly.

7.2.2. Cellular (Optional)

Cellular communication allows your therapy data to be wirelessly transmitted to your care provider to enhance the quality of your treatment. It also allows therapy settings to be updated in a more timely manner.

- (1). Make sure the Cellular function is on. (please refer to 4.4)
- (2). The icon displayed at the top right of your screen indicates the signal strength.
- (3). The data is automatically transmitted to the SleepWell Track after therapy has been completed. In order to make sure that your data is transferred, leave your device connected to the power outlet at all times and make sure that it is not in Airplane Mode.

Only for iX series 9S-010620

-  **NOTE: Once the Cellular mode is on, the device is automatically connected to the internet.**
-  **NOTE: Therapy data might not be transmitted if you use it outside of the country or region of purchase.**
-  **NOTE: Devices with cellular communication might not be available in all regions.**
-  **NOTE: Please be aware that within the wireless network, the availability and quality of the network may be affected by terrain, buildings, and the weather.**

8. Troubleshooting

The table below lists troubleshooting solutions for the problems that may occur. If the problem persists, contact your care provider's service agent or WELLELL. Do not open the device.

Problem	Possible Causes	Solutions
No display	<ul style="list-style-type: none"> ● The power cord is not connected to the power socket. ● LCD failure or controlled PCB failure. 	<ol style="list-style-type: none"> (1). Ensure the power cord is connected. (2). Contact your equipment provider for repair.
Illuminator under LCD is not on	<ul style="list-style-type: none"> ● Device is in sleep mode. ● LCD failure 	<ol style="list-style-type: none"> (1). Make sure the power cord is plugged in, then press any button/dial, and the LCD should automatically turn on. (2). If the above solution dose not work, contact your equipment provider for repair.
Buttons or Dial do not work	Button or Dial failure	Contact your equipment provider for repair.
Air delivery is slow	<ul style="list-style-type: none"> ● During ramp time. ● Filter is too dirty. ● Flow generator failure. 	<ol style="list-style-type: none"> (1). Check the ramp time setting (2). Change the filter regularly. (3). Contact your equipment provider for repair.
Data cannot be copied to the SD card	<ul style="list-style-type: none"> ● SD card is full. ● SD card is not inserted correctly. ● Data on the SD card is corrupted. 	<ol style="list-style-type: none"> (1). Ensure the SD card has enough capacity. (2). Ensure the SD card is inserted into the slot. (3). Format the card. If it dose not work, you might need to use a new SD card.
There is no air flow through the mask	<ul style="list-style-type: none"> ● Device is not turned on or working correctly ● Flexible tubing is not connected correctly. ● Air tubing is blocked. 	<ol style="list-style-type: none"> (1). Check the power cord is connected properly. (2). Reconnect the flexible tubing correctly. (3). Unblock the air tubing.
Condensation in mask or tubing	<ul style="list-style-type: none"> ● The Humidifier level setting is too high. ● The operating environment or position 	<ol style="list-style-type: none"> (1). Adjust to low temperature setting. (2). Remove any air conditioner nearby the heated

Problem	Possible Causes	Solutions
	of heated humidifier is not correct. ● The setting of heated tubing is not on.	humidifier. Or keep room temperature up near 25°C (3). Adjust up the level of heated tubing setting until the condensation reduced or become comfortable.
Water Leakage	● Water chamber is not properly assembled. ● Water chamber worn out. ● The air outlet of water chamber (humidifier) is not closed.	(1). Remove the water chamber from the heated humidifier, pour out the water and reassemble the water chamber again, making sure water chamber closes securely, and fill the water till it reaches maximum line and check if it still leaks or not. (2). Replace with a new water chamber. (3). Make sure Light grey tubing or heated tubing is applied on the air outlet of the water chamber (humidifier).
Air is leaking from around mask or tubing or device (Low Pressure) Screen: show flashing icon 	● Tubing is damaged. ● Tubing is not correctly connected. ● Humidifier is not properly assembled. ● Water chamber is not properly assembled. ● Mask is not properly connected. ● User does not wear the mask properly.	(1). Check the condition of the tubing. (2). Connect the tubing properly. (3). Connect the humidifier properly. (4). Connect the water chamber properly. (5). Make sure the mask is perfectly sealed on the user's face. (You can apply the "Mask Fit" function for checking your mask fitting condition.) (6). Screen: flashing  will disappear
Therapy data has not been sent to my care provider	● Wireless coverage may be poor. ● The No wireless connection icon is displayed on the top right of the screen. no	(1). Make sure that the device is placed where there is coverage (ie, on your bedside table, not in a drawer or on the floor). The icon of wireless should be

Problem	Possible Causes	Solutions
	wireless network available. ● Device may be in Airplane Mode. ● Data transfer is not enabled for your device.	turned on if the wireless coverage is good. (2). Turn off Airplane Mode (see 4.5) (3). Disconnect the device from mains supply, wait 5 seconds, then reconnect. The device will automatically re-upload data. (4). Make sure there are no co-existence interfere by adding new item/disconnect item to the connected IT-network. (5). Talk to your care provider about your settings

Caution:

- **Relevant risks need to be analyzed and evaluated by manufacturer. Relevant protective measures should be taken.**
- **Any changes to the relevant network configuration and network device connection methods may cause unknown risks. If necessary, please contact your care provider or manufacturer.**

8.1. Error and Alarm Messages shown on LCD

The table below lists Error and Alarm solutions for the issue that may happen. If the issue persists, contact your care provider’s service agent or WELLELL. Do not Open the device.

Message type	Message in LCD	Definition	Troubleshooting
Error Code	E001	Remove the Humidifier during treatment	Do not remove the humidifier during treatment.
	E002	SD communication fail	(1). Ensure the SD card is inserted into the slot. (2). Ensure the SD card has enough capacity. (3). Use SD card certified by Wellell.

Message type	Message in LCD	Definition	Troubleshooting
	E003	Remove the SD card while data is being processed	Do not remove SD while data is being processed.
	E004	Wi-Fi signal is too weak	<p>Please put the device closer to AP. Ensure the device is placed where there is coverage (ie, on your bedside table, not in a drawer or on the floor). The icon of wireless should be turned on if the wireless coverage is good.</p> 
	E005	No Wi-Fi connection	<p>(1). Ensure the Wi-Fi is on and Wi-Fi setting is correct. (or Turn off Airplane Mode, please see 4.5.)</p> <p>(2). Ensure the device is placed where there is coverage (ie, on your bedside table, not in a drawer or on the floor). The icon of wireless should be turned on if the wireless coverage is good.</p>  <p>(3). If the above solution is not working, please send the SD card to your care provider. The SD card also contains your therapy data.</p> <p>(4). Make sure there is no change in the IT-network configuration (or Make sure there are no co-existence interfere by adding new item/dis connect item to the connected IT-network)</p>

Message type	Message in LCD	Definition	Troubleshooting
	E006	Cellular is too weak	<p>Ensure the device is placed where there is coverage (ie, on your bedside table, not in a drawer or on the floor). The icon of wireless coverage should be turned on if the wireless coverage is good.</p> 
	E100-XX	System fault (primary function can be executed)	<ol style="list-style-type: none"> (1). Press any buttons or dial to return to standby mode. (2). Re-plug and the device will automatically re-upload data. (3). Contact your equipment provider. Do not open the device.
	E101-XX	System fault (primary function cannot be executed)	<p>Contact your equipment provider. Do not open the device.</p>
Alarm		Low Pressure	<ol style="list-style-type: none"> (1). Check the condition of the tubing. (2). Connect the tubing properly. (3). Connect the humidifier properly. (4). Connect the water chamber properly. (5). Close the chamber properly. (6). Make sure the mask is perfectly sealed on user's face.

9. Technical Specifications

9.1. Device

Item	Specifications
Power Supply	90W, DC 24V, 3.75A
Power Adapter	AcBel Polytech Inc., Model No.: ADH007 Input: 100-240V, 50/60Hz, 2.0-1.0 A Output rating: 24VDC, 3.75A, 90W
Mode of Operation	CPAP & APAP
Pressure Range	4~20 cmH ₂ O (adjustable in 0.2 cmH ₂ O increment)
Maximum Pressure	5~20 cmH ₂ O (adjustable in 0.2 cmH ₂ O increment)
Minimum Pressure	4~19 cmH ₂ O (adjustable in 0.2 cmH ₂ O increment)
Ramp Time	0~45 minutes (adjustable in 5-minute increment)
Ramp Starting Pressure	4 to therapy pressure cmH ₂ O (adjustable in 0.2 cmH ₂ O increment)
Operating Altitude	Up to 8,000 ft (2,438 m)
Altitude Compensation	Auto altitude compensation (up to 8,000 fts)
Communication	Wi-Fi 802.11b; 802.11g
	Bluetooth (for Wi-Fi configure setting) RF frequency: min 2379 MHz RF frequency: max 2496 MHz
	Cellular (optional) *Twelve Bands FDD-LTE: 700, 800, 850, 900, 1700/2100 (AWS), 1800, 1900, 2100, 2600 MHz (bands 1, 2, 3, 4, 5, 7, 8, 12, 18, 19, 20, 28) * Seven Bands UMTS (WCDMA/FDD): 800, 850, 900, 1700/2100 (AWS), 1800, 1900 and 2100 MHz (bands 1, 2, 4, 5, 8, 9, 19) *Quad Band GSM: 850, 900, 1800 and 1900 MHz Bd1, Bd2, Bd3, Bd4, Bd5, Bd7, Bd8, Bd12, Bd17, Bd18, Bd19, Bd20, Bd28

Item		Specifications												
		RF ID No <table border="1"> <tr> <td>9S-010600 / 9S-010640</td> <td>9S-010620 / 9S-010660</td> <td></td> </tr> <tr> <td>x</td> <td>FCC ID:QIPPLS62-W</td> <td>Cellular</td> </tr> <tr> <td>FCC ID:2AEQ402</td> <td>FCC ID:2AEQ402</td> <td>Wi-Fi</td> </tr> <tr> <td>FCC ID:2AEQ4RIFO</td> <td>FCC ID:2AEQ4RIFO</td> <td>BLE</td> </tr> </table> <p style="text-align: center;">  NOTE: The device should be used at a minimum distance of 8" (20 cm) from the body duration operation. </p>	9S-010600 / 9S-010640	9S-010620 / 9S-010660		x	FCC ID:QIPPLS62-W	Cellular	FCC ID:2AEQ402	FCC ID:2AEQ402	Wi-Fi	FCC ID:2AEQ4RIFO	FCC ID:2AEQ4RIFO	BLE
9S-010600 / 9S-010640	9S-010620 / 9S-010660													
x	FCC ID:QIPPLS62-W	Cellular												
FCC ID:2AEQ402	FCC ID:2AEQ402	Wi-Fi												
FCC ID:2AEQ4RIFO	FCC ID:2AEQ4RIFO	BLE												
Dimensions (W x D x H)		14.1 x 14.5 x 10.9cm; 23.5 x 14.5 x 14.7cm (with humidifier)												
Weight		Approximately 900 g; 1490g (with humidifier)												
Environment	Temperature	Operating: 5°C ~35°C (41°F ~95°F) Storage: -15°C ~50°C (5°F ~122°F) Transport: -15°C ~60°C (5°F ~140°F) <p style="text-align: center;">  NOTE: The device is immediately suitable for use if transported and stored according to the specified storage and transport conditions. </p>												
	Humidity	Operating: 15%RH ~ 95%RH non-condensing Storage: 10%RH ~ 90%RH non-condensing Transport: 10%RH ~ 90%RH non-condensing												
	Altitude	Operating: Sea level to 2,438 m; air pressure range 1060 hPa to 752 hPa												
Classification:		Class II Type BF, Applied Parts: Mask. <p style="text-align: center;">  Not suitable for use in the presence of a flammable anesthetic mixture. </p> IP22												
Air Tubing		WELLELL Light Gray Tubing – 15mm REF: SF00022 Flexible plastic, 1.83m												

Item	Specifications
	Operating/Storage/Transport conditions are the same as device.
Intake Filter (Ultra Filter)	WELLELL Material : Polypropylene (PP) Overall dimension: L 6mm x W 3mm Efficacy: NaCl Efficacy at 30 LPM: >93.66% (*Tested in accordance to TSI8130 NaCl .1 micron particle size) Pressure drop at 30 LPM: <0.7 mm H ₂ O (*Tested in accordance to TSI8130 pressure drop standard)
Service life	Device: 5 years Water Chamber: 2.5 years Light grey tubing: 1 year Heated tubing: 1 year
General	The patient is an intended operator.

NOTE: The manufacturer reserves the right to modify the specifications without notice.

9.2. Humidifier with Chamber for iX series (Optional)

Item	Specification
Model No.	9S-010800
Power Requirement	DC, 24V, 2.5A
Dimensions	11.5 x 14.2 x 14.7cm
Weight (without Water)	Approximately 590 g
Water Capacity	To maximum fill line 300ml
Maximum Heater Plate Temperature	Approximate 70°C (158°F)  NOTE: please do not touch the heated plate during/after operating.
Pressure Drop	0.2cmH2O @ 60LPM
Humidity Output	≥ 12 mgH2O/L (Ambient temperature: 23°C ± 2°C & Relative humidity: 60%±15%)
Operating Environment	Temperature: 5°C ~35°C (41°F ~95°F) Humidity: 15%-95%RH, non-condensing Altitude: Sea level to 2,438 m; air pressure range 1060 hPa to 752 hPa
Storage Environment	Temperature: -15°C ~50°C (5°F ~122°F) Humidity: 10%RH to 90%RH, non-condensing
Transport Environment	Temperature: -15°C ~60°C (5°F ~140°F) Humidity: 10%RH to 90%RH, non-condensing
Classification	Class II Type BF, Applied Parts: Mask  Not suitable for use in the presence of a flammable anesthetic mixture IP22

 **NOTE: The manufacturer reserves the right to modify the specifications without notice.**

9.3. Heated Tubing – 15mm for iX Auto series (Optional)

Item	Specification
REF No	SF00023
Internal Diameter:	15mm
Length	1.83 m
Temperature Range	1~6 levels, 18°C ~33°C (64.4°F ~91.4°F)
Cut-out	41°C (106°F)
Operating Environment	<p>Temperature: 5°C ~35°C (+41°F ~95°F) Humidity: 15 to 95% relative humidity, non-condensing Altitude: Sea level to 2,438 m; Sea level to 2,438 m; air pressure range 1060 hPa to 752 hPa</p> <p style="text-align: center;"> NOTE: The air flow for breathing produced by this therapy device can be higher than the temperature of the room. Under extreme ambient temperature conditions (104°F/40°C) the device remains safe.</p>
Storage Environment	<p>Temperature: -15°C ~+50°C (+5°F ~122°F) Humidity: 10 to 90% relative humidity, non-condensing</p>
Transportation Environment	<p>Temperature: -15°C to +60°C (+5°F to 140°F) Humidity: 10 to 90% relative humidity, non-condensing</p>
Resistance to Flow at rated flow of 60 l-min -1	≤ 0.1 kpa
Compliance Standards	<p>ISO5367:2014 ISO 80601-2-74:2017</p>

 **NOTE: The air tubing should be transported and stored in the original packaging and protected from dust and exposure to direct sunlight.**

 **NOTE: the manufacturer reserves the right to modify the specifications without notice.**

9.4. Noise Emissions Values

The A-weighted sound pressure level is:

Device: 28 dB(A) with and uncertainty of 2 dB(A).

Device with Humidifier: 29.5 dB(A) with and uncertainty of 2 dB(A).

The A-weighted sound power level is:

Device: 36.2 dB(A) with an uncertainty of 2 dB(A).

Device with Humidifier: 37.6 dB(A) with an uncertainty of 2 dB(A).

NOTE: Values determined according to noise test code given in ISO 80601-2-70:2015, using the basic standards ISO 3744 and ISO 7779.

9.5. Pressure Accuracy

Maximum static pressure variation at 10 cmH₂O (10 hPa) according to ISO 80601-2-70:2015

	15mm Light grey tubing W	15mm Heated Tubing W
Without Humidification	±0.5 cmH ₂ O (±0.5 hPa)	N/A
With Humidification	±0.5 cmH ₂ O (±0.5 hPa)	±0.5 cmH ₂ O (±0.5 hPa)

Static pressure accuracy has a measurement uncertainty of ± 0.15 cm H₂O

Maximum dynamic pressure variation according to ISO 80601-2-70:2015

Device without/with Humidification and 15mm Light grey tubing

Pressure(cmH ₂ O)	10 BPM	15 BPM	20 BPM
< 10 cmH ₂ O	0.65 / 0.65	0.8 / 0.8	1.6/ 1.8
≥ 10.0 to 20 cm H ₂ O	1/ 1.2	1.4 / 1.5	2.5 / 2.6

Device with Humidification and 15mm Heated Tubing

Pressure(cmH ₂ O)	10 BPM	15 BPM	20 BPM
< 10 cmH ₂ O	0.65	0.8	1.8
≥ 10.0 to 20 cm H ₂ O	1.2	1.5	2.6

Dynamic pressure accuracy has a measurement uncertainty of ± 0.98 cm H₂O

Results are expressed as STPD (Standard Temperature and Pressure, Dry).

9.6. Maximum Flow Rate

Pressure (cm H ₂ O)		Test pressures (cm H ₂ O)				
		4	8	12	16	20
15 mm tubing	Flow (LPM)	>100	>100	>100	>100	>100



Results are expressed as STPD (Standard Temperature and Pressure, Dry).

9.7. Maximum Single Fault pressure

When assessing the relative risks of using this equipment, the Physician should understand that this device can deliver pressures up to 20 cmH₂O. In the event of certain fault conditions, a maximum pressure is 40 cmH₂O. Results are expressed as STPD (Standard Temperature and Pressure, Dry).

10. Symbols

Symbol	Definition	Symbol	Definition
	BF symbol, which indicates this product is in accordance to the degree of protection against electric shock for type BF equipment.	IP22	Protected against solid foreign objects of 12.5 mm or greater and dripping water when tilted up to 15 degree.
	Follow Instructions for Use		Natural rubber latex is not present
	Refer to instruction manual	—MAX—	Maximum Fill Line (For Heated Humidifier)
	Caution In case of application as a safety sign, the rules according to ISO3864-1 are to be adhered to. See safety sign ISO7101-W001(Table D.2, safety sign 2)		Air flow direction
	Class II		Authorized representative in the European community

Symbol	Definition	Symbol	Definition
	Manufacturer		RTCA/DO-160 Section 21 Category M Compliant
	EC directive medical device label (Compliance to MDD 93/42/EEC; 4 digits number is NB's number)		Serial number
	Operating conditions		Storage conditions
	Hot Water, Hot Surface (For Heated Humidfier)		FCC logo
	JRF logo (for Japan only)		Stand-by
	Disposal of Electrical & Electronic Equipment (WEEE): This product should be handed over to an applicable collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.		

11. Servicing

The iX series systems are intended to provide safe and reliable operation according to the instructions provided by Wellell Medical. Wellell Medical recommends that the iX series systems be inspected and serviced by authorized technicians if there is any sign of wear or concern with device function. Otherwise, service and inspection of the devices generally should not be required during the five year design life of the device.

12. Limited Warranty

Wellell Inc. (here after 'Wellell') gives the customer a limited manufacturer warranty on new original Wellell products and any replacement part Autotod by Wellell in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. This warranty does not cover damage caused by accident, misuse, abuse, alteration and other defects not related to material or workmanship. To exercise your rights under this warranty, please contact your local, authorized Wellell dealer.

Product	Warranty period
CPAP devices (including external power supply units)	2 years
Heated humidifier	1 year
Water chamber	6 months
Mask systems	6 months
Light grey tubing (15mm)	6 months
Heated tubing (15mm)	6 months

Appendix A: EMC Information

Guidance and Manufacturer’s Declaration- Electromagnetic Emissions:

iX series is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment-Guidance
RF emissions CISPR 11	Group1	iX series uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	iX series is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	
<p>⚠ Warning:</p> <ol style="list-style-type: none"> iX series should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance (3.3m) between devices or turn off the mobile phone. 		

Guidance and Manufacturer’s Declaration- Electromagnetic Immunity:

iX series is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test Levels				
Basic EMC standard	Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT	Compliance Levels	Electromagnetic Environment-Guidance
Electrostatic Discharge (ESD) IEC61000-4-2	±8 kV contact ± 2, 4, 8, 15 kV air		±8 kV contact ± 2, 4, 8, 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	±2 kV, 100 kHz for power supply line ±1 kV, 100 kHz for input/output line		±2 kV, 100 kHz for power supply line ±1 kV, 100 kHz for input/output line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	± 0.5, 1 kV (line to line) ± 0.5, 1, 2 kV (line to earth)		± 0.5, 1 kV (line to line)	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	Voltage Dips: i) 100% reduction for 0.5 period, ii) 100% reduction for 1 period, iii) 30% reduction for 25/30 period, Voltage Interruptions: 100% reduction for 250/300 period		100 Vac & 240 Vac input	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	30 A/m	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in	Portable and mobile RF communications equipment should be used no closer to any part of

	<p>bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz</p>	<p>and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz</p>	<p>ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz</p>	<p>this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter.</p>
<p>Radiated RF EM Fields IEC61000-4-3</p>	<p>3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz</p> <p>380 - 5800 MHz, 9-28 V/m, Pulse modulation or FM ±5 kHz deviation , (1kHz sine)</p>	<p>10 V/m 80 MHz to 2,7 GHz 80 % AM at 1 kHz</p> <p>380 - 5800 MHz, 9-28 V/m, Pulse modulation or FM ±5 kHz deviation , (1kHz sine)</p>	<p>10V/m 80 MHz - 2.7 GHz 80 % AM at 1 kHz</p> <p>380 - 5800 MHz, 9-28 V/m, Pulse modulation or FM ±5 kHz deviation, (1kHz sine)</p>	<p>Recommended separation distance</p> <p>$d = \sqrt{P}$ 150kHz to 80MHz $d = 0.6\sqrt{P}$ 80MHz to 800MHz $d = 1.2\sqrt{P}$ 800 MHz to 2.7GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).^b</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey , should be less than the compliance level in each frequency ranged.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

<p>NOTE 1: U_T is the a.c. mains voltage prior to the application of the test level</p> <p>NOTE 2: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people</p> <p>NOTE 4: The performance of wireless functions in the presence of Electromagnetic Interference has not been evaluated.</p>
<p>a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.</p> <p>b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.</p>

Recommended separation distances between portable and mobile RF communications equipment and this device:

iX series is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \sqrt{P}$	80 MHz to 800 MHz $d = 0.6\sqrt{P}$	800 MHz to 2.7 GHz $d = 1.2\sqrt{P}$
0.01	0.1	0.06	0.12
0.1	0.31	0.19	0.38
1	1	0.6	1.2
10	3.1	1.9	3.8
100	10	6	12

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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All rights reserved
Wellell America Corp.