

INNOVATIVE S-SERIES SCOOTERS

From the shopping mall to the country, the new luxury S-Series by STERLING are ready for any journey. Packed full of brilliant features with superb driving dynamics, and looking sharp from every angle, there's an S-Series mobility scooter to suit everyone.













INNOVATION IN **DESIGN**

With advanced LED lighting technology, the S-Series scooters are 400 times more energy efficient than standard bulbs, meaning more battery power and range every trip.



→ Innovative active anti-tip wheels

Move up as you ride over obstacles,
and keep you safe on slopes.

Safety as standard, S-Series models all feature rear high visibility LED indicators, position and brake lights.



SZZOOCOMPACT SCOOTER

The S400 **6kph** luxury compact scooter with big features, is highly maneuverable around shopping centres or busy high streets, and with the all-around suspension as standard makes every journey a comfortable one.

With **38 Ah** batteries the S400 is ideal for trips around town, plus with the bright, energy efficient LED's lighting your way, you're just as visible on those late-night shopping trips as you are during the day.











Max Speed: **6 kph**



Colour: Silver Metallic



Batteries: 38 Ah



Wheel size: 25 cm



Max. user weight: 136 kg



CONTROLS DESIGNED FOR YOU

Finding your perfect driving position is a breeze with the infinitely adjustable tiller. Clever design means the headlight automatically tilts as the tiller is adjusted.



S425 COMPACT OUTDOOR SCOOTER

The luxury S425 crams in all of the features of its little brother, but with added power!

Taking the bright, front/rear LED lights, comfortable all-around suspension and combining it with a maximum speed of **12 kph**, the S425 is suitable for road use as well as on pavements.

The increased wheel height will help to tackle those usual outdoor obstacles, and longer journeys can be handled with ease as the larger **55 Ah** batteries will help to power your way.











Max Speed: 12 kph



Colour: Carbon Metallic



Batteries: **55 Ah**



Wheel size: **30 cm**



Max. user weight: **150 kg**







COMFORT ON EVERY TRIP

The S-SERIES advanced seating solution, takes comfort to the highest level, from the seat height, depth and recline adjustments, to the flip up, width, angle and depth adjustable comfort armrests — every feature has been designed for comfort — for you.



→ Infinite seat depth adjustable



→ Quick backrest angle adjustment



→ Seat rotation enables comfortable transfers

S/OO PERFORMANCE OUTDOOR SCOOTERS

The luxury S700, with superior outdoor performance is your perfect companion for outdoor adventures, to take you where you want to go.

There's no need to worry about distance either as the large **75 Ah** battery size will power an extended range and a lively top speed of **15 kph**.

The high performance design means gradients, and higher kerbs are tackled with ease and the large tyres and advanced 4-wheel suspension provides a smooth, rewarding drive each and every time.











Max Speed: 6 kph / 12 kph / 15 kph



Colour: Silver-/Carbon-/
Black Metallic



Batteries: **75 Ah**



Wheel size: **33 cm**



Max. user weight: 160 kg

TECHNICAL DATA

	N177-30	The state of the s		
S-SERIES	<i>5</i> 400	<i>5</i> 425	5700	
Max speed	6 kph	12 kph	6 kph / 12 kph / 15 kph	
Max. user weight	136 kg	150 kg	160 kg	
Battery size	38 Ah x 2	55 Ah x 2	75 Ah x 2	
Charger	5 A	8 A	8 A	
Max. potential range*	36 km	42 km	45 km	
Max. power (peak)	1000 W	1200 W	2700 W	Ca Carlot
All around LED lighting	Yes	Yes	Yes	
Easy adjust tiller & wig wag control	Yes	Yes	Yes	
Suspension	Yes	Yes	Yes - advanced	Parket
Tyres	25 cm pnuematic	30 cm pnuematic	33 cm pnuematic	Special date of
Overall width / length	64 cm / 134 cm	66 cm / 135 cm	67 cm / 150 cm	attern)
Seat width x depth	46 cm x 46 cm	49 cm x 49 cm	49 cm x 49 cm	
Seat height range	45 cm - 50 cm	45 cm - 50 cm	45 cm - 50 cm	
Minimum turning distance	1.9 m	2.0 m	2.2 m	
Colour	Silver Metallic	Carbon Metallic	Silver / Black / Carbon Metallic	
Scooter weight inc. batteries	109 kg	117 kg	146 kg	

^{*} Max potential range based on industry benchmark range figures for similar products and battery sizes

All over Europe, thousands of people have made their lives easier with a STERLING mobility scooter by Sunrise Medical. For over 30 years, STERLING have designed and built mobility scooters to the highest standards, giving the comfort and safety every customer needs. For further information on product options please check the order form on our website. You can view this document also in a larger font on a PDF-file at www.SunriseMedical.eu







Sunrise Medical GmbH & Co. KG

Kahlbachring 2-4 69254 Malsch/HD

Tel.: +49 (0) 72 53 / 9 80-0 Fax: +49 (0) 72 53 / 9 80-222

E-Mail: international@sunrisemedical.de

www.SunriseMedical.eu