


TERUMO NEEDLE AGANI[™]



**A wide range of needles for beautiful performance
from 30G to 18G**

TERUMO NEEDLE AGANI 30G - 18G features:

- High-quality, stainless-steel structure assures smooth puncture.
- Large internal diameter delivers high-volume flow.
- Semi-transparent needle hub is convenient for blood reflux observation.
- Meticulously designed, sharpened needle point ensures smooth puncture, reduces tissue damage and minimizes pain.
- Needle hub color coding conforms to international gauge size standards, for clear and instant recognition of needle size.

TERUMO NEEDLE AGANI product range

Gauge	Color	Length	Bevel	Product Code
18G	pink	38 mm (1 1/2")	Regular	AN*1838R
19G	cream			AN*1938R
20G	yellow	25 mm (1")		AN*2025R
		38 mm (1 1/2")		AN*2038R
21G	green	25 mm (1")		AN*2125R
		38 mm (1 1/2")		AN*2138R
22G	black	25 mm (1")		AN*2225R
		38 mm (1 1/2")		AN*2238R
23G	blue	25 mm (1")		AN*2325R
		32 mm (1 1/4")		AN*2332R
		38 mm (1 1/2")		AN*2338R
24G	purple	25 mm (1")		AN*2425R
		38 mm (1 1/2")		AN*2438R
25G	orange	16 mm (5/8")		AN*2516R
		25 mm (1")		AN*2525R
26G	brown	13 mm (1/2")		AN*2613R
27G	grey		AN*2713R	
30G	light-yellow		AN*3013R	



MANUFACTURED BY : Zhejiang Kindly Medical Devices Co.,Ltd

NO.252, Yongqiang Road, Yongzhong, Longwan, Wenzhou, Zhejiang Province, P.R.C.

SOLD BY : TERUMO CORPORATION

44-1, 2-chome, Hatagaya, Shibuya-ku, TOKYO, 151-0072, JAPAN

© Registered Trademark
AGANI is a trademark of TERUMO CORPORATION
© TERUMO CORPORATION, March 2012
SS30E-312KL-IO(3.12)A